

Nestin Ab, chicken polyclonal, Novus Biologicals, NB100-1604

Used on mouse neural stem cells (2 confocal images). Cells fixed with 4% PFA.

Blocking: 5% NGS, 0.1% Triton in PBS, 30 min.

Ab incubation: 1/1000 dilution in 1% NGS, 0.05% Triton in PBS, O/N at 4C.

Detection: anti-chicken alexa IgG-546 fluorescent antibody (1/1000 dilution).

Incubated in 1% NGS, 0.05% Triton in PBS, 2h at RT.