


PKN3 Antibody Characterization

Comments:

- Total protein in the listed cell lines or transfected HEK293 cells was analysed using the Bradford assay and protein samples (20µg each) were resolved by SDS-PAGE. Immunoblots were screened:
 - *anti*-PKN3 (Rabbit, NBP1-30102, Novus) antibody (1:4000 dilution)
- NBP1-30102 detecting FLAG-PKN3 transfected HEK293 cells. This antibody is entirely specific to PKN3. The expression of these proteins confirmed using anti-FLAG. This antibody also readily detects the respective endogenous protein.



*Research provided by Dr. Eamon Mulvaney, University College Dublin, Ireland